

## P-915K Low-Profile Piezo Objective Scanner For High Scanning Frequencies



- Very Low Profile of 15 mm
- Travel Range 75  $\mu\text{m}$
- Clear Aperture for Objectives with W0.8 x 1/36" Thread
- Frictionless, High-Precision Flexure Guiding System for Better Focus Stability and Minimized Runout
- Very Low Profile
- Outstanding Lifetime Due to PICMA<sup>®</sup> Piezo Actuators

Model	Active axes	Travel range	Resonant frequency @ 150 g	Dimensions
P-915KLPZ Objective Scanner	Z	75 $\mu\text{m}$	200 Hz	60 x 60 x 15 mm